
## Deutsch

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
## English

### 1. Important information - read prior to start-up

You have acquired a valuable Riester product manufactured in compliance with Directive 93/42/EEC for medical products and subject to continuous stringent quality control. Should you have any queries, please contact the Company or your Riester Agent who will be pleased to assist you. For addresses see last page of these Operating Instructions. The address of your authorised Riester Agent will be supplied to you on request. Please note that any products described in these Operating Instructions are only suited for application by trained operators. Please also note that correct and safe operation of products will only be guaranteed when Riester products and accessories are used throughout.

**Warning:** Never connect the pneumatic tourniquet to pressure-generating equipment.

#### 2. Purpose

The Riester pneumatic tourniquet was designed for restricting the blood circulation to limbs in case of amputations.

#### 3. Components (see figure)

- Manometer
- Pump
- Cuff
- Release valve
- Tube connection part I (male)
- Tube connection part II (female)

#### 4. Operation of the unit

- Close the release valve (3.4.) by turning it clockwise
- Screw tube connector at the end of the cuff tube into the tube socket on the manometer housing.
- The pneumatic tourniquet with cuff must be checked before each application by trained staff or a physician on function and density.

#### 5. Applying the cuff

Place the cuff around the respective limb and tighten the straps on the cuff by pulling them through the metal buckles. The tube exit should be approx. in the centre of the limb to be treated.

#### 6. Replacement of bladder

Remove the bladder from the side opening of the cuff and insert the new bladder through the opening. Pull the tube through the opening on the long side of the cuff.

#### 7. Operating the pump

Hold the pump (3.2.) in one hand and, using the other hand, carefully pull the air pipe of the pump on the black plastic ball (3.9.) up to the stop and then push it back into the valve. Repeat the pumping action until the manometer scale shows the desired value.

#### 8. Restricting the blood circulation to the limbs

- Restricting the blood circulation to the arms In order to restrict the blood circulation to the arms, the manometer must be pumped up to approx. 250 mmHg, using the pump.
- Restricting the blood circulation to the legs In order to restrict the blood circulation to the legs, the manometer must be pumped up to approx. 550 mm Hg with the pump.
- Pressure control The adjusted pressure which the manometer shows must be continuously supervised by a responsible physician during the application and, if necessary, be corrected by after-pumps.

#### 9. Letting out the air from the cuff after use

Open the deflating valve (3.4.). The air can now escape from the cuff and the cuff can be removed.

#### 10. Cleaning Cleaning and/or disinfecting General information

The goal of cleaning and disinfection of medical products is the protection of patients, users and third persons and conserving the value of the medical products. On account of the product design and the used material, no defined limit of maximum processing cycles can be fixed. The lifetime of the medical products depends on their function and on a appropriate treatment of the devices. Before returning faulty products for repair they must have gone through the described reprocessing process.

**Cleaning:** Manometer can be cleaned outside with a humid cloth until optical cleanliness is achieved.

**ATTENTION!** Never place the manometer in liquid! This item is not approved for automated reprocessing and sterilization. These procedures cause irreparable damage!

**Cleaning:** After removing the bladder, wipe the cuff covers with a damp cloth. Alternatively, these can be washed with soap and cold water. If you decide on the latter course, rinse the cuffs with clear water afterwards and let them air dry. Wipe the bladder and tubes with a damp cloth.

**Disinfection:** After removing the bladder, wash the cuff covers in cold water to which disinfectant has been added. Afterwards, let them air dry. Only disinfectants with approved efficiency and in accordance with the national state dards can be used. The bladder and tubes can be wiped with a cotton cloth moistened with ethanol.

**IMPORTANT!** Do not iron the cuff covers. Never expose the cuffs to intensive solar radiation! Never touch the cuff covers or bladders with a sharp instrument, since this could cause damage!


#### 11. Spare parts

No. 11221 Manometer  
No. 11220 Pump  
No. 11222 Upper arm cuff with bladder  
No. 11227 Cuff cover for upper arm bladder  
No. 11224 Pad for upper arm cuff  
No. 11219 Children's cuff with bladder  
No. 11229 Cuff cover for children's cuff  
No. 11246 Pad for children's cuff  
No. 11223 Thigh cuff with bladder  
No. 11228 Cover for thigh cuff  
No. 11225 Pad for thigh cuff

#### 12. Maintenance/accuracy test

The pneumatic tourniquet and its accessories do not require any special maintenance.

For the accuracy test, remove the tube from the manometer and hold the manometer in vertical position. If the pointer points to the zero position on the scale, the unit has been set up correctly. If the pointer does not point to the zero position, please return the unit to ourselves or your local authorised Riester agent, whose details are available from ourselves.

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## WARRANTY

This product has been manufactured under the strictest quality standards and has undergone a thorough final quality check before leaving our factory. We are therefore pleased to be able to provide a warranty of

**2 years from the date of purchase** on all defects, which can verifiably be shown to be due to material or manufacturing faults. A warranty claim does not apply in the case of improper handling. All defective parts of the product will be replaced or repaired free of charge within the warranty period. This does not apply to wearing parts. For R1 shock-proof, we grant 5 years, for precisa N shock-proof, 3 years for the calibration, which is required by CE-certification. A warranty claim can only be granted if this Warranty Card has been completed and stamped by the dealer and is enclosed with the product. Please remember that all warranty claims have to be made during the warranty period. We will, of course, be pleased to carry out checks or repairs after expiry of the warranty period at a charge. You are also welcome to request a provisional cost estimate from us free of charge. In case of a warranty claim or repair, please return the Riester product with the completed Warranty Card to the following address:

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